



**BUREAU
VERITAS**

MINERALS

▶ SS80, SS230, SSXX

Package Description	Soil Preparation
Samples Digestion	NA
Instrumentation Method	NA
Legacy Code	SS80, S230, SSXX
Applicability	Soils and Sediments

▶ METHOD DESCRIPTION

Wet or damp soil samples are dried at 60°C (Air dried or 40°C if specified by the client). Soil and sediment sieved to -80 mesh (SS80) or -230 mesh (SS230), unless client specifies otherwise (SSXXX). Sieves cleaned by brush and compressed air between samples.

NOMINAL OPENING	US STANDARD SIEVE	TYLER SIEVE
250 µm	60	60
180 µm	80	80
150 µm	100	100
106 µm	140	150
75 µm	200	200
63 µm	230	250

